



STANLEY

Handling your safety



Jolly

mobility stair climber

The Jolly is a tracked mobile stair climber which allows the user to travel up and down the stairs in the comfort of their own wheelchair. Designed to manage the majority of staircases, it is compatible with electric chairs, wheelchairs with tilting or movable backrests and even pushchairs.

The unique design keeps the user at a level position on the stairs, instilling confidence and trust at all times during operation.

The light frame can easily be dismantled into two separate pieces to be loaded and transported by car.


Permanent installation is not required as the wheelchair is simply attached to the tracked base and a trained operator may then simply manoeuvre the user up and down the stairs in a smooth motion.

Benefits

- Tracked mobile stair climber provides a smooth motion during the climbing process
- Easy dismantle for vehicle loading
- Loading ramps provided
- Accommodates various wheelchair types
- Travel in the comfort of your own wheelchair
- Bespoke attachments available
- Overcomes architectural challenges
- Practical, manoeuvrable and light
- Simple push button control
- Free demonstrations
- National maintenance coverage
- Certificated training provided
- Competitive finance and hire options
- 24-month parts warranty included

Specification

- Maximum capacity: 150Kg
- Adjustable headrest and safety belt
- Key switch for authorised personnel use
- Emergency button
- Non-skid and non-marking tracks
- Double security device
- One person emergency descend function
- Indicator for checking stairway angle
- Ascent range 50 floors with a weight of 150Kg
- Motor: 12V – 160W, direct current
- Battery: 12V 24Ah, sealed, rechargeable, low maintenance



Jolly ramp



STANLEY HANDLING

Unit 3, Sovereign Park,
Laporte Way, Luton,
Bedfordshire, LU4 8EL

t: 0800 298 2980

e: sales@stanleyhandling.co.uk

w: www.stanleyhandling.co.uk

