



STANLEY

Handling your safety



NuDrive Air

lever drive technology

NuDrive is the innovative lever-drive propulsion accessory for manual wheelchairs, providing more mobility whilst improving physical well-being. With NuDrive you don't need to propel your chair by gripping the wheel rims by hand; instead, you can propel yourself forwards, backwards, manoeuvre and brake, simply by pushing the levers.

The 'Snap-Lock' technology allows each of the wheel adaptor levers to be fitted in seconds, without altering your current wheel configurations. Once in place, the lightweight design adds no width to your chair and can be left unobtrusively positioned on the wheel. When you want to use NuDrive, the drive units can simply be clicked into place and then removed and placed back into the carrying case when finished.

Benefits

- Reduces the force needed to self-propel by up to 50%
- Protects hands and wrists, and improves hygiene
- Propulsion enforces a healthier lifestyle
- Reduces the risk of shoulder degradation and injury by shifting shoulder loading
- Overcome elevations and inclines with greater ease
- Easily attaches to most manual wheelchair types with 22, 24 and 25-inch spoked wheels
- Great for outdoor use on normal & difficult terrain
- Snap-lock levers are fitted in seconds
- Levers are counter-balanced meaning they stay in an upright position at all times
- Secure control of the wheels, even when wet
- Safer control on declines
- Wrist and hand pain reduced in over 65% of users
- Improved breathing due to better posture
- Extends the range and ability of your manual wheelchair

Specification

Aluminium wheel adaptor

- Added width: 0cm
- Weight: 1Kg each
- Colour: metallic grey powder coated

Aluminium lever drive units

- Added width: 1.5cm each, assuming 3 degree wheel camber
- Weight: 2.5Kg each
- Max user weight 136Kg



STANLEY HANDLING

Unit 3, Sovereign Park,
Laporte Way, Luton,
Bedfordshire, LU4 8EL

t: 0800 298 2980

e: sales@stanleyhandling.co.uk

w: www.stanleyhandling.co.uk

